



Inside Out Technique for Central Venous Occlusions

Monnie Wasse, MD, MPH

Vice Chair & Professor of Medicine

Director, Interventional Nephrology

Rush University Medical Center

May 11, 2019

VASA Practicum, Houston, TX

COI


- Principle investigator for Surfacer U.S. IDE clinical study

CVO is a Significant Problem in Vascular Access



An estimated **2M people** ww require treatment with dialysis or kidney transplant to stay alive¹

Permanent arteriovenous fistulas (AVFs) remain the preferred form of vascular access




Yet, **central venous catheters (CVCs)** serve as an **important bridge** while waiting to secure access

⚠ Chronic CVC use leads to central venous obstruction (CVO) in 40% of cases²

¹<https://www.kidney.org/kidneydisease/global-facts-about-kidney-disease>

²Data on file at Bluegrass Vascular Technologies, Inc.



BVT has pioneered an innovative inside-out approach to achieve reliable and repeatable RIJ access for hemodialysis patients.

Surfacar System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

Technique video available at <http://bluegrassvascular.com/surfacar/>

No Current Standard of Care

- No standard care algorithm when confronting occluded central veins
 - Re-establish access in the current vein or “try another vein”
 - Over time, all central veins become destroyed
- Lack of solution drives limited appreciation of problem and downstream consequences of progressively occluding veins

The Surfacar® Story

A UNIQUE APPROACH TO RESTORE VASCULAR ACCESS

Developed by **interventional cardiologist John Gurley, M.D.** at a leading U.S. academic institution

Over **600 Inside-Out®** procedures performed to-date

CE Mark Approved
July 2016

Dr. John Gurley, M.D., Interventional Cardiologist
Founder Of Bluegrass Vascular Technologies




RESTORE ACCESS.
PRESERVE OPTIONS.

Surfacer® Inside-Out® Access Catheter System

Inside-Out® Approach

- Unique, minimally invasive approach
- Delivered by inserting a guidewire, under fluoroscopy, through the femoral vein




Surfacer System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

Surfacer® Inside-Out® Access Catheter System

Targeted Navigation

- The Surfacer® is navigated up through the vasculature of the torso
- The device then exits out a desired exit point in the RIJ



Surfacer System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

Europe

CE Mark received 3Q'16

Commercial Feedback Forms

(Field Evaluations)

27 hospitals
90+ patients

SAVE Registry

Post-Market Clinical
Follow-up

5 sites
30 patients

United States

Under investigational use only


SAVE-US IDE Study

7 sites
30 patients

when completed, clinical data on >150 cases






*Surfacer System is CE Marked – Not available for sale in the United States. Limited by
Federal (or United States) law to investigational use.*

Surfacer Procedure Time



Surfacer System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

SAVE Registry Occlusion Types*



- Bilateral IJ and SC occlusions (Type 1)
- Brachiocephalic (Type 2)
- SVC occlusion above azygos (Type 3)
- Total occlusion of the SVC (Type 4)

Surfacer System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

SAVE Registry Results

- All endpoints met (primary and secondary)
 - Primary Performance (successfully place access): 97%
 - Primary Safety: (acute safety, absence of device related AEs): 100%
 - Secondary Performance (ability to place access): 97%
 - Secondary Safety (technique conversion): 97%
- No adverse events or post-op complications reported
- Fluoroscopy times averaged 6.8 min
- Contrast use on average: 29.7 cc
- Average procedure times
 - Workstation Sheath Insertion to Removal = 23 min (6-70 min range)
 - Surfacar Insertion to Needle Wire Exit = 5 min (0-60 min)

Surfacar System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.


Surfacer System US IDE Clinical Study

- IDE study Inclusion: up to 30 patients that require central venous access
- Primary endpoints:
 - ▶ Performance:
 - Rate of successful transient CVCs created across venous occlusions where success definition is placement of patent operational CVCs.
 - ▶ Safety:
 - Acute device safety, defined as the absence of procedural complications at discharge and 7 days post-procedure
 - Overall device and procedure related anticipated adverse events compared to historical safety data from CV placement procedures.
- Secondary endpoints:
 - ▶ Performance:
 - Device able to advance from femoral vein to subclavian exit to facilitate catheter placement, etc.
 - ▶ Safety:
 - Technique conversion rate/causes and catheter malposition

Surfacer System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.

IDE Surfacr System Interim Results

- **28** patients enrolled, **25** treated successfully with 2 conversions due to tortuous anatomy and 1 technique aborted due to abnormal anatomy
- 100% access gained/device performed in treated patients, with no serious device related adverse events
- All primary and secondary endpoints met



Surfacr System is CE Marked – Not available for sale in the United States. Limited by Federal (or United States) law to investigational use.